
METRIA

RSD205 series

Please read the User Manual carefully before use, and follow all operating and safety instructions!

user manual

english / español / français

User Manual

RSD Portable Digital Refractometer

Preface

Users should read this manual carefully, follow the instructions and procedures, and beware of all the cautions when using this instrument for a maximum performance and elongate the equipment life.

Service

When help is needed, you can always contact the Service Department of manufacturer for technical support:

www.labbox.com / e-mail: info@labbox.com

Please provide the customer attention service representative with the following information:

- Serial number
- Description of problem
- Your contact information

Warranty

This instrument is warranted to be free from defects in materials and workmanship under normal use and service, for a period of 12 months from the date of invoice. The warranty is extended only to the original purchaser. It shall not apply to any product or parts which have been damaged on account of improper installation, improper connections, misuse, accident or abnormal conditions of operation.

For claim under the warranty please contact your supplier.

1. Safety Instructions

**Warning:**

- Read the operating instructions carefully before use
- Ensure that only trained staff works with the instrument

- If the equipment is used in any manner not specified by the manufacturer, the protection provided by the equipment may be affected.
- We are not responsible for any injury or damage caused by modifying this instrument without authorization, nor for any injury or damage from the use of the instrument for any purpose not foreseen.

2. Inspection

2.1 Reception inspection

Unpack the equipment carefully and check for any damages which may have arisen during transport. Please contact your supplier for technical support.

**Note:**

If there is any apparent damage to the system, contact your dealer.


2.2 List of items

The package includes the following items:

Items	Qty
Refractometer	1
Batteries	2
Carrying case	1
Soft cleaning cloth	1
User Manual	1

Table 1

3. Control elements

4. Control panel

Item description	Use
 POWER	Press to turn on the device or turn it off by holding button for 3 seconds
 ZERO	Press to set calibration at zero or hold 2 seconds to switch between °C and °F.
 READ	Press to read the samples or hold it during 2 seconds to switch between different measuring scales.
 POWER + READ	Reset to factory setting

5. Calibration / Zero point calibration

MODEL	SCALES	RANG	DIVISIONS	ACCURACY	CALIBRATION
RSD200	Brix	0-45%	0.1	±0.1	Deionized water
	Refractive index	1.3330-1.4098	0.0001	±0.0003	
RSD500	Brix	0-85%	0.1	±0.5	Deionized water
	Refractive index	1,3330-1,5100	0.0001	±0.0005	
RSD205	Oe Germany		1	±1	Deionized water
	Oe (Swiss)	3-150 / 0-150	1	±1	
	Mass S/W	0-35% / 0-22%	0,1	±0,2	
	Vol AP	0-25	0,1	±0,2	
	KMW		0,1	±0,2	

The zero position on the refractometer calibration must be checked before each use and afterwards. A zero point calibration should be performed before any measurement where the highest possible accuracy is required or if the environmental temperature has great variations. It is recommended to use deionized water. (The temperature for the calibration liquid should be ideally near 20°C / 68°F)

- Inspect the prism to ensure that is clean and dry.
- Place a few drops of calibration liquid on the prism

- Press **ZERO** (calibration button) to start
- The refractometer will show the calibration process on the LCD (Image 1)
- If the procedure is successful a "PASS" message will appear on the screen (Image 2)
- After the calibration, the LCD will return to temperature mode. Calibration result will be saved and it will be the new zero point after the device is turned off and on again.

Note:

- After calibration, please remember to clean the liquid from the prism surface and keep the prism surface dry and clean to avoid erosion.


Image 1


Image 2

6. Measurement

- Make sure that the prism surface is clean and dry.
- Place a few drops of the sample on the prism.
-  Press, **READ** and “---” will appear on the screen.
- The result(s) will be stored in the device for up to 60 seconds. (Image 3)
- Press the “Power” button to re-check the previous measurements.
- After testing, please remember to clean and dry the surface to maintain the device in optimal condition

Image 3

7. Scale Selection

-  Hold **READ** during 2 seconds, and the scale will change (Image 4)
- Repeat the step 1 until the desired scale is shown on the screen
- The refractometer will save the last selected scale.

Image 4

8. Temperature selection

-  Hold **ZERO** during 2 seconds, and the temperature unit will switch between Celsius (°C) and Fahrenheit (°F)
- The refractometer will save the last selected temperature unit.

9. Resetting

- Press  +  and the option for returning to factory settings will be shown on the screen (Image 5)
- Press  to confirm resetting, or press  to cancel.

Image 5

10. Power on/off

1. Press  to turn on the device.
2. Press  to turn off the device, LCD will show the process (Image 6)
3. To be energy efficient, the device will enter into sleep mode after 60 seconds of inactivity. And will be active after the power button is pressed again.
4. The refractometer will turn off after 90 seconds of no operation.

Image 6

11. Trouble shooting

Zero setting error (Image 7)

- Out of calibration range.

Action:

- Make sure you are using the correct calibration solution, referred on the specification table.
- Make sure there is sufficient amount of calibration liquid to cover the prism surface.


Image 7

Measurement Error

- Exceeds the specified measurement range (Image 8)
- The reading is below the specified measurement range (Image 9)
- Not enough liquid on the prism (Image 10)

Image 8

Image 9

Image 10

Temperature error

- Temperature exceed range (Image 11)

Action:

- Check temperatures of sample liquid, device and ambience.

- Temperature under range (Image 12)

Action:

- Check temperatures of sample liquid, device and ambience.

Image 11

Image 12

Battery error

- Low battery indication (Image 13)

Action:

-Replace the batteries when the icon is red. 2 x AAA batteries are needed.

Image 13

12. Maintenance

Please keep the main surface clean and dry. It is extremely important to clean it after each use with a wet, soft, clean cloth or paper, dampened with a mild liquid soap and water. This may prevent cross contamination between samples and provide accurate subsequent reading. The use of solvents or petroleum-based cleaners is not recommended.

Important Precautions:

- To ensure that the LCD displays properly, do not expose the instrument to a low or high temperature or prolonged exposure to direct sunlight.
- Prevent the equipment from getting wet and suffering blows that could affect the equipment during transport.
- Do not disassemble or assemble the instrument or change the inner parts.
- To avoid accuracy being affected by evaporation, ensure that the measurement is realized immediately after the solution has been dropped into the device.
- Do not use the instrument under a high humidity or corrosive environment.
- Avoid strong light sources during the measurements (such as sunlight, lamp etc.).
- When storing the instrument for long periods of time, it is recommended to remove the batteries. Use only AAA batteries. Pay close attention to battery polarity, reversing the polarity can damage the instrument.

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