

Technical Specification

Heating Mantle Body	Made of Metal (Aluminum Alloy/MS depending in size)
Heating Mantle Material	Glass Fibre Yarn
Mantle Capacity	To hold standard Round Bottom Flasks upto 20000ml
Temperature Range	Approx. 400°C
Heating Control	Energy power regulator
Permissible Ambient Temp.	5- 40 °C
Permissible relative moisture	80%
IP Class	10

DAN
— LOGITECH —

Heating Mantle

OPERATING MANUAL

1000.DNEU

1000.DNUS



Dan Logitech Instruments Pvt. Ltd.

#282, HSIIDC Ind. Estate, Saha, Ambala Cantt - 133104, Haryana (INDIA)

info@danlogitech.com | www.danlogitech.com

Congratulations on your purchase of a DAN Logitech Heating Mantle.

General Description:-

Heating Mantle -The heater are made of fiber glass yarn with heating elements covered in fiber Glass sleeve, to provide uniform heating of Borosilicate flask up to 400°C.The mantle sizes are as per round bottom DIN 12347/ISO 4797 Standard.

Safety Information:- This product has been designed for safe operation when used as detailed in accordance with the Manufacturer's instructions.

General Remarks:-

- The operator has to make sure that the Heating mantle is stable.
- Put only Glass flask on the heating mantle that fit the heating mantle's volume.
- Always replace the fuse from the same type & current.
- Always put a Glass flask in the heating mantle. Never put beaker & anything else or never operate without flasks.
- Never rotate the Heat knob out of their ranges to avoid damages.
- These Heating mantle should not be operated unattended.
- Do not lift unit by the top.
- Don't spill any liquids into the heating mantle.

Important Safety Instruction:-

- Never turn on the Heat knob without flask.
- The Glass flask should match the size of the heating mantle exactly, it needs good heat conductor between the vessel & mantle.
- Heating mantle are to be operated in accordance with individual national regulations or the specifications of the professional associations in each case.
- Keep the heating mantle dry. Even small quantities of liquid can destroy the heating elements.
- These heating mantles are not moisture protected, and must therefore be protected against humidity and damp, in the particular against aggressive media a few drops of which could be sufficient to destroy the heat conductor. A heating mantles which has become wet may only be used again once it has been thoroughly dried out for half an hour (for preference in a drying cabinet).
- For safety reasons heating mantle should be operated with appropriate temperature limiters.
- Do not touch the heating surface when they are switched on.
- Do not touch the surface of the heating element while energized or hot.
- The upper area of the heating mantle can also take high temperature and must therefore not be touched during operation. Do not use metal instruments near


Model No. (230 VAC/50Hz)	Model No. (115 VAC/60Hz)	Mantle Capacity (ml)	Wattage (W)
1000.DNEU.01	1000.DNUS.01	50	150
1000.DNEU.02	1000.DNUS.02	100	150
1000.DNEU.03	1000.DNEU.03	250	150
1000.DNEU.04	1000.DNUS.04	500	200
1000.DNEU.05	1000.DNUS.05	1000	300
1000.DNEU.06	1000.DNUS.06	2000	450
1000.DNEU.07	1000.DNUS.07	3000	450
1000.DNEU.08	1000.DNUS.08	4000	450
1000.DNEU.09	1000.DNUS.09	5000	600
1000.DNEU.10	1000.DNUS.10	6000	600

Maintenance & Repair:-

According the safety regulations the operator has to verify regularly the safety and functioning of the heating mantle (just as other electronic devices, repair must be done at the factory or by authorized companies.)

Operating Instructions :-

Over view of heating mantle heat control



Correct Usage :-

1. Place the unit on the level table or bench near to the mains supply
2. Put the flask with liquid inside the heating mantle. Make sure that the flask makes good contact with the yarn.
3. Connect the instrument to a ground mains outlet, after ensuring that the Voltage is the same as given on sticker. The unit must be connected to a single phase A.C power supply with protective earth.
4. Ensure that the heat knob is at OFF position.
5. Connect the heating mantle to power source.
6. Set the required temperature with the heat knob by rotating.
7. LED indication glows which shows that the heating has started.
8. After use keep the unit in a safe place.
9. When the heating is complete, switch the heat by knob to their off positions. Disconnect the mains electricity supply.

Features:-

1. Mantle knitted with glass fiber yarn.
2. Provides uniform heating to flasks.
3. Heats up to 400°C.
4. For round bottomed Flasks as per DIN-12347/ISO 4797.
5. Inbuilt protection for liquid spill.
6. Easy maintenance.
7. Long life.

the heating element surface of the mantle.

- Do not install the heating mantle in hazardous area. They are not explosion protected and accordingly may not be used for heating & stir processes in which there may be a risk due to explosive media or explosive gas air mixtures.
- To avoid electrical shock, completely cut off power to the unit by disconnecting the power cord from the unit or unplug from the wall outlet. Disconnect unit from the power supply prior to maintenance and servicing.
- To avoid electric shock always use a properly grounded electrical outlet of correct voltage and current handling capacity.
- You can use a soft moistened cotton cloth to clean the outer surface of the heating mantle.
- Don't use in the presence of pressurized or sealed containers. Fire or explosion may result, causing severe injury or death.

Symbols On The Instruments:-

⚠ Caution:-for continued safety replace only with fuse of same rating.

⚠ Caution:- Hot surface

Unpacking & Contents:-

Please unpack the equipments carefully and check for any damages. It is important that any damages which may have arisen during transport are ascertained when unpacking. If applicable a fact report must be sent immediately.

- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Socket must be earthed (protective ground contact).
- Attention - Magnetism! Effects of the magnetic field must be taken into account (e.g. data cardiac, carriers paceman keys...).
- Risk of burns! The unit can reach temperatures in excess of 400°C. Pay attention the residual heat after switching off. Please make sure that the mains cable does not touch the hot surfaces.

- Wear your personal protective equipment in accordance with the hazard category of the medium to be processed. Otherwise there is a risk of: splashing liquids - projectile parts release any toxic or combustible gases.
- Set up the appliance in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- The feet of the appliance must be clean and undamaged.
- Position the knob at the left stop before starting up. Gradually increase the speed.
- Reduce the speed if - the medium splashes out of the vessel because the speed is too high-the appliance is not running smoothly - the container moves on the set-up surface.
- Caution! Be aware of the flash points of media to be used power must always be set to at least 25 °C lower than the fire point of the media used.
- When using PTFE-coated magnetic bars, the following has to be noted:
Chemical reactions of PTFE occur in contact with molten or dissolved alkaline and alkaline - earth metals, as well as with fine-particle powders of metals of the 2. and 3. group of the periodical system at temperatures above 300-400°C. Only elementary fluorine, chlorine trifluoride and alkaline metals do attack PTFE, halogen hydrocarbons have a reversibly swelling effect.
- Check the appliance and accessories beforehand for damage each time you use them. Do not use damaged components.
- Only replace damaged parts with spare parts identical to the original in function and quality.
- Beware of the risk of Contents- flammable materials-glass breakage as a result of mechanical shaking power-incorrect container size-too much medium unsafe condition of container
- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced

in other ways, e.g. through light irradiation.

- Do not operate the appliance in explosive atmospheres with hazardous substances or under water.
- A separation from the line is made with the equipment only by pulling net and/or device plug.
- Safe operation is only guaranteed with the accessories described in the “Accessories” chapter.
- Always disconnect the plug before fitting accessories.

Please check the contents of your carton against the relevant product .

Item No.	Description	Qty
1.	Heating mantle	1
2.	User manual	1

Installation:-

Electrical Safety And Installation

This equipment is designed to be used safely under the following conditions:-

- Indoor use.
- Mains supply voltage fluctuations up to $\pm 10\%$ of the nominal voltage. This equipment must be used with proper earthed / grounded to a fixed earth / grounded mains socket outlet. The mains supply is to earthed / grounded in accordance with current legislation. See Technical Specification for recommended fuse ratings.
- Ensure only the correct rated mains input fuses are fitted. (Where applicable ensure the correct Mains cord and moulded IEC plug and lead set fuse if fitted).
- Check the voltage on the product data label on this product unit and those of any accompanying electrical accessory. Ensure the rating conforms to your local supply.
- Do not install this product or accessories on a surface which may become flooded.

Observation: - The surface of the heating element of a mantle cartridge will upon receipt look slightly discoloured. This discolouration is normal and occurs at the factory during test when the mantle is first heated up.