
H03D Digital magnetic stirrer with heating
S03 Magnetic stirrer without heating
S03D Digital magnetic stirrer without heating

Please read the User Manual carefully before use, and follow all operating and safety instructions!

user manual

english / español / français

User Manual

H03D Digital magnetic stirrer with heating
S03 Magnetic stirrer without heating
S03D Digital magnetic stirrer without heating

Preface

Users should read this Manual carefully, follow the instructions and procedures, and beware of all the cautions when using this instrument.

Service

If help is needed, you can always contact your dealer or Labbox via www.labbox.com

Please, provide the customer service representative with the following information:

- Serial number
- Description of the problem
- Your contact information

Warranty

This instrument is guaranteed to be free from defects in materials and workmanship under normal use and service, for a period of 24 months from the date of invoice. The warranty is extended only to the original purchaser. It shall not apply to any product or parts which have been damaged on account of improper installation, improper connections, misuse, accident or abnormal conditions of operation.

For claim under the warranty, please contact your supplier.

1. Safety Instructions

	Warning! <ul style="list-style-type: none">• Read the operating instructions carefully before use.• Ensure that only trained staff works with the instrument.
	Risk of burn! <ul style="list-style-type: none">• Caution when touch the housing parts and the hotplate, which can reach temperatures of 280°C.• Pay attention to the residual heat after switching off the equipment.
	Protective ground contact! <ul style="list-style-type: none">• Make sure that the socket is earthed (with protective ground contact) before use.

- When working, wear personal protective equipment to avoid the risks from:
 - Splashing and evaporation of liquids
 - Release of toxic or combustible gases
- Set up the instrument on a spacious, stable, clean, non-slip, dry, and fireproof surface. Do not operate the instrument in explosive atmospheres, with hazardous substances, or under water.
- Gradually increase the speed. Reduce the speed if:
 - Stirring bar breaks due to high speed
 - The instrument is not running smoothly or the container moves on the base plate
- Temperature must always be set at least 50 °C below the flash point of the media used.
- Be aware of hazards due to:
 - Flammable materials or media with a low boiling temperature
 - Overfilling of media
 - Unsafe container
- If the stirring bar contains PTFE, please note:
 - Elemental fluorine, trifluoride ion, and alkali metals will corrode PTFE. Halogen hydrocarbons expand PTFE at room temperature.
 - Molten alkali, alkaline earth metals or their solution, as well as elements on second and third column of the Periodic Table of elements will react with PTFE when temperature reaches 300 ~ 400 °C.
- Process pathogenic materials only in closed vessels.
- Check the instrument and the accessories prior to each use. Do not use damaged components. Safe operation is only guaranteed with the accessories described in the “Accessories” chapter. Accessories must be securely attached to the device and cannot come off by themselves. Always disconnect the plug before the assembly of any accessory.
- If the external temperature sensor is needed, the tip of the measuring sensor must be at least 5-10 mm away from the vessel bottom and wall.
- The instrument can only be disconnected from the main power supply by pulling out the mains plug or the connector plug.
- The voltage stated on the label must correspond to the main power supply.
- Ensure that the cable does not touch the hotplate. Do not cover the device.
- The instrument may only be opened by experts.

- Keep away from high magnetic fields.
- Observe the minimum distances between devices, between the device and the wall, and above the assembly (min. 100 mm)


Figure 1

2. Proper Use

The instrument is designed for mixing and/or heating liquids in schools, laboratories, or factories. This device is not suitable for use in residential areas or other restrictions mentioned in Chapter 1. **Not tested for continuous work (24h/7days).**

3. Inspection

3.1 Unpacking

Unpack the equipment carefully and check for any damages which may have arisen during transportation. Please contact the manufacturer/supplier for technical support.

Note:

If there is any apparent damage to the system, please do not plug it into the power line.

3.2 Listing of Items

The package includes the following items:

Items	Qty
Main unit	1
Power cable	1
Stirrer bar	1
User Manual	1

Table 1

4. Control

4.1 Control elements




Figure 2: H03D series

Speed Control Knob


LED Display

Figure 3: S03D series

Speed Control Knob

Figure 4: S03 series

Power switch

Figure 5: S03, S03D, H03D series

Connection PT1000

lbox Socket
instruments

Figure 6 : H03D series

	Items	Descriptions
H03D series	Temperature Control Knob "Heat"	Set the temperature parameters. The function "heating" is switched ON or OFF by pushing the ON/OFF knob. To increase the value turn the knob clockwise.
	Speed Control Knob "Stir"	Set the rotary speed. The function "Stirring" is switched ON or OFF by pushing the ON/OFF knob. To increase the value turn the knob clockwise.
	LED Display	The LED screen displays the real working state and all the settings.
	Probe	When the external temperature sensor PT1000 is plugged in, the sensor icon  lights up.
	Power Switch	Switch the instrument ON or OFF.
S03D series	Temperature Control Knob "Stir"	Set the rotary speed. The function "Stirring" is switched ON or OFF via push ON/OFF knob. To increase the value turn the knob clockwise.
	LED Display	LED displays the real working state and all settings.
	Power Switch	Switch the instrument ON or OFF.
S03 series	Speed Control Knob "Stir"	The stirring function is switched ON by rotating the knob. To increase the value turn the knob clockwise. The stirring function is switched OFF by rotating the knob all the way clockwise.
	Power Switch	Switch the instrument ON or OFF.

Table 2


Figure 7: H03D series

Characters	Descriptions
Temperature display area	When heating function is switched ON, LED displays the temperature setting value and shifts to real value in 5 seconds. When the heating function is switched OFF and the hotplate temperature is still above 50°C, LED displays “Hot”, otherwise LED displays “OFF”.
Speed display area	When stirring function is switched ON, LED displays the speed setting value and flashes. The setting value does not flash until real speed reaches the set point.

Table 3


Figure 8: S03D series

Characters	Descriptions
Speed display area	While stirring function is on, the LED screen displays the speed setting value and flashes. The set value does not flash until the real speed reaches the set point.

Table 4

5. Trial Run

- Make sure the required operating voltage and power supply voltage match.
- Ensure that the socket is properly grounded.
- Plug in the power cable, ensure the power is on, and begin initializing.
- Add the medium into the vessel with an appropriate stirring bar.
- Place vessel on the work plate.
- Set the target stirring speed and begin.
- Observe the stirring bar and LED display (S03D series and H03D series).
- Set the target temperature and start heating.
- Observe the LED display (H03D series).
- Stop the heating and stirring functions.

If the operations above cause no abnormalities, the

have been damaged during transportation, please contact manufacturer/supplier for technical support.

Warning!

Forbidden to transfer the vessel while the instrument is working.

6. Function: Heating (H03D series)

The device is controlled by digital temperature control technology, which has two separate safe circuits. The hotplate is kept at a constant temperature by a digital control circuit. The hotplate temperature can also be monitored from a separate, adjustable safe circuit. The two internal temperature sensors for temperature control are built into the hotplate. The single external PT1000 can monitor the temperature of the sample.

- If using an external sensor PT1000: Always connect it to the device BEFORE turning the power on and insert it into the sample BEFORE heating (see section 6.1).
- Set the temperature by rotating the temperature control knob slowly towards the target value.
- While the heating function is on, the LED screen displays the temperature value on the left-hand side.
- The heating function is switched on or off by pushing the heating knob.

The heating function will be stopped automatically under abnormal conditions. The safe circuit keeps the hotplate temperature under 320°C.

Please, in that case, follow the instructions below:

- Switch OFF the instrument.
- Ensure the external temperature sensor is inserted in the media heated.
- Switch ON the instrument and run heating function.

If the heating function does not work, please contact manufacturer/supplier for technical support.

Once turned on, the instrument automatically displays the last running speed and temperature parameters.

Generally, the LED screen cannot display the actual temperature of sample in the vessel or hotplate surface. Temperature differences exist between:

- Hotplate center and outer edge.
- The sample in the vessel and hotplate surface.

In order to ensure the accuracy of the temperature inside the container, please use the external temperature sensor PT1000.

6.1 Working with external temperature sensor


Figure 9

The external temperature sensor PT1000 is the manufacture's standard accessory. Compared to the temperature control of the hotplate, the external temperature sensor can control the medium's temperature more precisely.

- **If an external temperature sensor is needed, it must be plugged in before turning on the device** (fig. 6).
- Ensure the external temperature sensor is connected firmly and inserted into the media before heating.
- Switch on the instrument and start heating.
- The LED screen displays the temperature setting value for 5 seconds and then shifts to the real value.

Warning!

Forbidden to plug in/pull out external temperature sensor in heating.

6.2 Residual heat warning (HOT)

In order to prevent burns from the hotplate, the digital hotplate has a residual heat warning function. When the heating function is switched off and the heating plate temperature is still above 50°C, "Hot" will flash on the LED screen to warn that there is a hazard of burns from the hotplate. When the unit is powered off, the LED screen alternates the temperature of the hotplate and "Hot". When the hotplate temperature drops to below 50°C, the unit will automatically switch off. If users want to turn off the LED immediately, just pull out the plug directly. When the plug is pulled out, the residual heat warning function stops running.

7. Function: Stirring

The function "stirring" of LED digital model (H03D and S03D):

- Turn on the device by pressing the switch ON/OFF.
- Start the function "stirring" by pressing the "STIR" knob and adjust the speed slowly until the set value.
- When starting the function "stirring", the display shows set value and then blinks. This blinking will stop only when the actual speed reaches the set value.

- At the end of the manipulation, turn the knob to the left to slow down and turn off "stirring" by pressing the knob.
- Turn off the device by pressing the switch ON/OFF.

The function "stirring" of classic model (S03):

- Turn on the device by pressing the switch ON/OFF.
- Activate and adjust the speed of stirring by slowly turning the speed control to the set value.
- At the end of the manipulation, turn the knob fully counterclockwise to turn off the stirring.
- Turn off the device by pressing the switch ON/OFF.

8. Faults

- Instruments cannot be power ON.
 - Check whether the power line is unplugged.
- Fault in power on self-test
 - Switch OFF the unit and then switch back ON.
- Stir speed cannot reach set point
 - Excessive medium viscosity may cause abnormal speed reduction of the motor
- Unit cannot be powered off by pressing the switch on / off.
 - Check if the residual heat warning function is still ON and hotplate temperature is above 50°C (the LED screen still works and "Hot" appears on the screen).

If these faults are not resolved, please contact manufacturer/supplier.

9. Maintenance and Cleaning

- Proper maintenance can keep instruments working properly and lengthen their lifetime.
- Do not spray cleanser onto the instrument when cleaning.
- Unplug the power line when cleaning.
- Only use recommended cleansers:

Dyes	Isopropyl alcohol
Construction materials	Water containing surfactants / Isopropyl alcohol
Cosmetics	Water containing surfactants / Isopropyl alcohol
Foodstuffs	Water containing surfactants
Fuels	Water containing surfactants

Table 5

- Before using any other cleaning or decontamination method, the user must verify with the manufacturer that the method will not harm the instrument. Wear proper protective gloves during the cleaning of the instrument.

Note:

- Electronic devices cannot be cleaned with cleanser.
- If maintenance service is needed, the instrument must be cleaned in advance to avoid pollution from hazardous substances, and sent back in the original packing.
- If the instrument will not be used for a long time, please switch off and place in a dry, clean location, at room temperature, and on a stable surface.

10. Standards and Regulations

Construction in accordance with the following safety standards:

EN 61010-1
UL 3101-1
CAN/CSA C22.2(1010-1)
EN 61010-2-10

Construction in accordance with the following EMC standards:

EN 61326-1

Associated EU guidelines:

EMC-guidelines: 2004/108/CE
Instrument guidelines: 73/023/EWG

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference. In that case, the user will be required to correct the interference at his own expense.

11. Specifications

Items	Specifications	
	H03D / S03D series	S03 series
Voltage [VAC]	*100~120/200~240 100~240	
Frequency [Hz]	50/60	50/60
Power [W]	*515/15	10
Stirring point position quantity	1	1
Max. stirring quantity (H ₂ O) [l]	3	3
Recommended magnetic bar [L, mm]	30 - 50	30 - 50
Motor type	DC motor	
Max. power input of motor [W]	5	5
Max. power output of motor [W]	3	3
Speed range [rpm]	100 – 1500, increment:10	0 – 1500
Rotary speed display	LED	Scale
Plate material	*Ceramic coated/Plastic	
Dimensions of workplate (mm)	Ø135	Ø135
*Heating power [W]	500	-
*Temperature range [°C]	RT-280, increment:1	-
*Temperature display [°C]	LED	-
*Temperature display accuracy [°C]	1	-
*The safe temperature of hotplate [°C]	320	-
*Temperature sensor in medium	PT1000	-
*Control accuracy of heating temperature with temperature sensor [°C]	±0.5	-
*Residual heat warning	50°C	-
Dimensions (mm)	*220×160×95 220×160×75	
Weight [kg]	*1.4/0.7	0.7
Permitted ambient temperature [°C]	5 - 40	
Permitted relative humidity	80%	
Protection class acc. to DIN 60529	*IP21/IP42	

*H03D series

Table 6

12. Accessories

Please contact our company to order the following accessories available for these devices: magnetic stirring bars, PT1000 temperature sensor with and without glass coated, and fixation support for the temperature sensor.

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