
## **H01 Magnetic stirrer with heating**

*Please read the User Manual carefully before use, and follow all operating and safety instructions!*

# **user manual**

english / español / français / italiano

# User Manual

## H01 Magnetic stirrer with heating

### Preface

Users should read this Manual carefully, follow the instructions and procedures, and beware of all the cautions when using this instrument.

### Service

If help is needed, you can always contact your dealer or Labbox via [www.labbox.com](http://www.labbox.com)

Please, provide the customer service representative with the following information:

- Serial number (on the back side)
- Description of the problem
- Your contact information

### Warranty

This instrument is guaranteed to be free from defects in materials and workmanship under normal use and service, for a period of 24 months from the date of invoice. The warranty is extended only to the original purchaser. It shall not apply to any product or parts which have been damaged on account of improper installation, improper connections, misuse, accident or abnormal conditions of operation.

For claim under the warranty, please contact your supplier.

# 1. Safety Instructions

	Connect the device to a grounded power supply to ensure the safety of the machine and the experiment. Connect the power as the machine requires it.
	The use of this equipment in inflammable and explosive, poisonous, or highly corrosive experiments is forbidden.
	Put the magnetic stirrer on a horizontal, flat, stable table.
	This equipment must only be used by previously qualified staff.
	When heating flammable substances, the hot plate temperature must be lower than the substance's flash point.
	Act with the utmost care when working with liquids that reach temperatures higher than 100°C (for example, oil).
	If the vessel is overfilled, it may cause overheating of some working parts, which will dissolve the flammable liquids.
	While the machine is working, do not touch the hot plate, nor the vessel or the solution in order to prevent high-temperature burns.
	Read the instruction book before operating the equipment.

Table 1

- When using the equipment, wear personal protective equipment to avoid the risks from:
  - Splashing and evaporation of liquids
  - Release of toxic or combustible gases.
- Set up the instrument on a spacious, stable, clean, non-slip, dry, and fireproof surface. Do not operate the equipment in explosive atmospheres, with hazardous substances or under water.
- Gradually increase the speed. Reduce the speed if:
  - The stirring bar breaks due to high speed
  - The equipment is not running smoothly or the container moves on the base plate.
- Temperature must always be set at least 25°C below the flash point of the media used.
- Beware of hazards due to:
  - Flammable material or media with a low boiling temperature
  - Overfilling of the container
  - Unsafe container
- Pathogenic materials must always be processed in close vessels.

- If the stirring bar contains PTFE, please note:
  - Elemental fluorine, trifluoride ion, and alkali metals will corrode PTFE. Halogen hydrocarbons expand PTFE at room temperature.
  - Molten alkali, alkaline earth metals or their solution, as well as elements on second and third column of the Periodic Table of elements will react with PTFE if the temperature reaches 300 ~ 400 °C
- Check the instrument and the accessories for damage prior to every use. Do not use damaged components. A safe operation is only guaranteed with the accessories described in the “Accessories” chapter. Accessories must be securely attached to the device and cannot come off by themselves. Always disconnect the plug before the assembly or disassembly of accessories.
- Ensure that the external temperature sensor is inserted into the media to a depth of at least 20 mm.
- When using metal vessels, do not place the temperature sensors on the bottom of the vessel. Placing sensors on the vessel bottom can cause temperatures that are too high to be measured, especially in media that have poor conductivity. The tip of the measuring sensor must be at least 5 mm away from the vessel bottom; a distance of 10 mm is ideal.
- The instrument can only be disconnected from the main power supply by pulling out the mains plug or the connector plug.
- The voltage stated on the label must correspond to the main power supply.
- Ensure that the cable does not touch the heating base plate. Do not cover the device.
- Keep away from high magnetic fields.
- Observe the minimum distances between devices, between the device and the wall, and above the assembly (min.100mm).


Figure 1

## 2. Proper Use

The instrument is designed for mixing and/or heating liquids in schools, laboratories, or factories. This device is not suitable for use in residential areas or other restrictions mentioned in Chapter 1.

### 3. Inspection

#### 3.1. Receiving Inspection

Unpack the equipment carefully and check for any damages that may have arisen during transportation. Please contact the manufacturer/supplier for technical support.

**Note:**  
If there is any apparent damage to the system, please do not plug it into the power line.

#### 3.2. Items list

The package includes the following items:

Items	Qty
Main unit	1
Power cable	1
Stirrer bar	1
User Manual	1
Rack with rod	1

Table 2

### 4. Description and Control

Figure 2: H01 series


Figure 3: H01 series

## 5. Trial Run

- Make sure the required operating voltage and power supply voltage match.
- Ensure that the socket is properly grounded.
- Plug the power cable into the system, then into the outlet. Make sure the power is on and begin the initiation.
- Add the medium into the vessel with a magnetic stirring bar; make sure the magnetic stirring bar has a convenient size.
- Put the vessel on the work plate.
- Set the stirring speed and start stirring.
- Observe the magnetic stirring bar.
- Set the target temperature and start heating.
- Observe the real temperature through an external control of the medium's temperature.
- Stop the heating and the stirring.

*If the operations above cause no abnormalities, the device is ready to operate. Otherwise, the device may have been damaged during transportation, please contact the manufacturer/supplier for technical support.*

## 6. Operation

- Place the magnetic stirrer on top of a horizontal, flat, and stable table making a free space of at least 30 cm on each side. Do not put the stirrer near any warm supply.
- Plug in the power cable and press the general switch.
- In the middle of the plate, place a glass with the magnet and the solution to be stirred; then, gradually increase the stirring speed.
- Setting the stirring speed: Turn the speed knob clockwise; the movement indicator will light up. Adjust the stirring speed according to your needs.
- Setting the temperature: Turn the temperature knob clockwise; the temperature indicator will light up. Then, adjust temperature's intensity according to your needs. Remember that the solution's temperature is always lower than the hot plate's adjusted temperature.
- Placing the thermometer: Assemble the two parts of the stirring bar and insert it in its support at the back side of the equipment. Once the bar is placed, place the bosshead at the desired height, hold it with the screw, and place the thermometer clamp.

**Note:**

- Increase the stirring speed gradually.  
Reduce the speed if:
  - Too high speed makes the stirrer break.
  - The equipment’s functioning is irregular or it is not standing stable on the heating plate.
- Make sure the power supply and voltage are enough for the system.
- Make sure the power line is at a safe distance from the heating plate.
- Do not cover the equipment.
- In case the equipment breaks down, cut off the power supply right away.

## 7. Faults

- Instruments cannot be switched ON
  - Check whether the power line is unplugged.
- Stirring speed cannot reach set point
  - Excessive medium viscosity may cause abnormal motor speed reduction.
- Heating does not after setting the temperature or stirring does not start by adjusting the control knob.
  - Check the control panel has damages which may have arisen during transport.

*If these faults are not resolved, please contact the manufacturer/supplier.*

## 8. Maintenance and Cleaning

- Proper maintenance can keep instruments working properly and lengthen the lifetime.
- Do not spray cleanser onto the instrument when cleaning.
- Unplug the power line when cleaning.
- Only use recommended cleansers:

Dyes	Isopropyl alcohol
Construction materials	Water containing surfactants / Isopropyl alcohol
Cosmetics	Water containing surfactants / Isopropyl alcohol
Food products	Water containing surfactants
Fuels	Water containing surfactants

Table 3

- Before using any other cleaning or decontamination methods, the user must verify with the manufacturer that the method will not harm the equipment. Wear proper protective gloves while cleaning the equipment.
- If maintenance service is needed, the equipment must be cleaned to avoid contamination from hazardous substances and shipped in its original packing carton.
- Use the instrument in a dry, clean, environment with a stable temperature.

## 9. Storage and transportation

- Keep it in a dry and clean room with good ventilation and no corrosive gas.
- Prevent it from wetting by the rain and avoid violent collision during transportation.

## 10. Associated Standards and Regulations

Construction in accordance with the following safety standards:

EN 61010-1  
UL 3101-1  
CAN/CSA C22.2(1010-1)  
EN 61010-2-10

Construction in accordance with the following EMC standards:

EN 61326-1

Associated EU guidelines:

EMC-guidelines: 89/336/EWG  
Instrument guidelines: 73/023/EWG

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

**NOTE:** This equipment has been tested to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

## 11. Specifications

Model	Voltage (V)	Speed (r.p.m.)	Heating Power (W)	Max. Vol. (ml)	Max. Temp. (plate)	Hot Plate size (mm)
H01	220V/50Hz	100 – 2000	180	2000	380	120 x 120

Table 4

## 12. Working condition

Ambient temperature: 5~40°C;

Ambient humidity: ≤90%;

Voltage: 220V±10%, 50/60Hz or 110V+/-10%, 50/60Hz

## 13. Accessories

Please contact our company to order the following accessories available for these devices: magnetic stirring bars.

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