## WB01 Thermostatic Water Bath

*Please read the User Manual carefully before use, and follow all operating and safety instructions!*

# user manual

english / español / français

# User manual

## WB01 Digital thermostatic water bath

### **Preface**

Users should read this Manual carefully, follow the instructions and procedures, and beware of all the cautions when using this instrument.

### **Service**

If help is needed, you can always contact your dealer or Labbox via [www.labbox.com](http://www.labbox.com)

Please, provide the customer service representative with the following information:

- Serial number (on the back side)
- Description of the problem
- Your contact information

### **Warranty**

This instrument is guaranteed to be free from defects in materials and workmanship under normal use and service, for a period of 24 months from the date of invoice. The warranty is extended only to the original purchaser. It shall not apply to any product or parts which have been damaged on account of improper installation, improper connections, misuse, accident or abnormal conditions of operation.

For claim under the warranty, please contact your supplier.

# 1. Safety Instructions

	Connect the device to an earthed power supply to ensure safety of the machine and the experiment; connect the power as the machine requires it.
	The use of this instrument in inflammable, explosive, poisonous, or highly corrosive experiments is forbidden.
	Place the water bath on a horizontal, flat, stable table leaving 30 free cm on each side.
	The heating mantle must only be used by previously qualified staff that has read the instructions manual and knows how to operate it.
	Do not place the heating mantle near any heat source.
	During its functioning, dangerous materials such as flammable or pathological substances must be out of the device's safety area.
	An overfilling of the vessel could cause some working parts to overheat, which could dissolve flammable materials and start a fire.
	While the machine is working, do not touch the working surface, in order to avoid burns.
	Read the instructions manual before using this device.

- When working, wear the necessary personal protective equipment to avoid the risk of:
  - Burns caused by splashing and evaporation of liquids
  - Intoxication caused by release of toxic or flammable gases.
- Set up the instrument on a spacious, stable, clean, non-slip, dry, and fireproof surface that can support the weight of the equipment. Do not operate the instrument in explosive atmospheres or with hazardous substances.
- The equipment must be used only when the resistor is covered by water; doing it otherwise will harm the resistor.
- The temperature must be set at least 25°C below the flash point of the substances used.
- Beware of hazards due to:
  - Flammable materials or media with a low boiling temperature
  - Overfilling of the vessel

- Unsafe vessels
- Process pathogenic materials only in closed vessels.
- Check the instrument and accessories before handling for damage prior to every use. Do not use damaged components. Safe operation is only guaranteed with the accessories described in the “Accessories” chapter. Accessories must be securely attached to the device and cannot come off by themselves. Always disconnect the plug before the assembly or disassembly of accessories.
- Use original spares only, consult your supplier for more information.
- The instrument can only be disconnected from the main power supply by pulling from the plug, not the cable.
- The voltage stated on the label must correspond to the main power supply.
- Ensure that the mains cable does not touch the heating base. Do not cover the device.
- Keep away from high magnetic fields.
- Pay attention to the temperature setting. Never leave the equipment unattended while the heating function is ON.

## 2. Using instructions

The instrument is designed for heating liquids in schools, laboratories or factories. It is not suitable for domestic use or for use in environments that can be hazardous for either the user or the instrument.

## 3. Inspection

### 3.1 Unpacking

Unpack the equipment carefully and check for any damages that may have arisen during transportation. If necessary, contact your supplier for technical support.

**Note:**

If there is any apparent damage to the equipment, please do not plug it into the power line.

### 3.1. Items list

The package includes the following items:

WB01 series:

Items	Quantity
Main unit	1
Instructions manual	1
Power cable	1
Resistor covering platform	1
Lid with hole	1
Fuse	1
Reducing discs	Según Modelo

## 4. Control Elements

4.1. Thermostatic water bath, WB01 series:

In the TWBA-001-001 series, the control panel is set up the following way:

1. "Power": ON/OFF switch
2. Control panel
3. Power inlet

In the other models, the control panel is set up this way:


- 1.1 'Power': ON/OFF switch
- 1.2 Drain switch: 'Outlet'
- 1.3 Control panel
- 1.4 Draining pipe
- 1.5 Power inlet

### Control Panel

The control panel is the same for all models.

- 2.1 SET: parameters confirmation
- 2.2 ◀: shift digits
- 2.3 ▼: decrease parameter
- 2.4 ▲: increase parameter
- 2.5 PV screen: measurement value
- 2.6 SV screen: parameter value
- 2.7 RUN: work indicator
- 2.8 OUT: heat indicator
- 2.9 AT: auto setting indicator
- 2.10 ST: time setting indicator
- 2.11 ALM: alarm indicator

## 5. Trial run

- Place the equipment on a flat, stable surface that can support the weight of the device leaving 30 free cm on each side.
- Make sure the required operating voltage and the power supply's voltage match.
- Make sure the socket is earthed.
- Place the resistor covering platform and fill the water barrel 2/3 of its total capacity.
- The use of decalcified water is recommended.
- Connect the power cable to the equipment and then to the socket. Make sure the power is ON.
- Turn ON the equipment by pressing the ON/OFF switch.
- The PV screen will display the real temperature of the water. If this is on, the electronic equipment is working properly.
- The SV screen will display the previously set temperature. Press the SET button and set a higher temperature than the water's real temperature. If the temperature increases, the equipment is working properly.

*Note: DO NOT use purified or demineralized water since these are great oxidants.*

## 6. Operating

*Note: before using the water bath, check the water level inside the barrel; if necessary, add more water. Remember that running the resistor without being completely covered may damage it.*

- Place the equipment on a flat, stable surface that can support the weight of the device leaving 30 free cm on each side.
- Press the green ON/OFF switch at the front part of the water bath. In the TWBA-001-001 series, the switch is on the left side.
- Once the equipment is ON, you will hear a light sound; this indicates that the equipment is ready to use.
- The PV screen displays the real water temperature, while the SV screen displays the previously set temperature.
- After using, turn off the equipment by pressing the ON/OFF switch and unplug it. Remember that, if the equipment is not going to be used for a long time, it must be drained. This function is not available in the TWBA-001-001 series; therefore, you need to wait for the equipment and the water to cool down before draining the water bath.

### 6.1. Set the working temperature

The max working temperature is 100°C; set the temperature this way:

- Press SET. The PV screen will display the command "Sy" and the SV screen will display the temperature value to be adjusted.
- Use the ▲ and ▼ buttons to set the desired temperature. Use the ◀ button to shift digits.
- Press SET to confirm the command; the resistor will start working to reach the set temperature. The PV screen will display the real temperature while the SV screen will display the set temperature.
- In order to get a homogeneous temperature in all the water, stir regularly. Adapt the reducing discs to the container's size.

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## 6.2. Timer setting

By default, every time the water bath is switched ON, this feature will be OFF. Therefore, the equipment will keep working until it is manually switched OFF by the user with the ON/OFF switch.

- Press the SET button for 3 seconds. The PV screen will display the “Loc” command and the SV screen will display, by default, the digits “0000”.
- Using the ▲ and ▼ buttons, configure the “0002” setting and press SET. The PV screen will display the command “St” and the SV screen will display “Off”.
- Use the ▲ and ▼ buttons (and the ◀ button to shift) to configure the time period in which you want to work. The maximum time value is 9999 min.
- Press SET to save the time value, and then wait to go back to the initial screen where the real temperature is displayed on the PV screen and the time, on the SV screen.

When the timer option is on, the indicator RUN will stay lit (green) and the ST and OUT indicators will flash (green); and once the time is over the RUN and OUT indicators will turn off and the ST indicator will keep flashing.

Once you turn ON this option, it will stay ON. To turn it OFF, proceed the following way:

- Press the SET button for 3 seconds; the PV screen will display the “Loc” command and the SV screen will display, by default, the digits “0000”. Then, configure the option “0002” and press SET. To deactivate the timer, set the time to “0000”; “Off” will appear on the SV screen. Press SET to save and wait to go back to the initial display.

*Note: before programming the timer, you must adjust the temperature you want to work with.*

## 6.3. Alarm setting

If the temperature of the water bath is 10°C above the set temperature, an alarm will sound, while the PV screen displays the real temperature of the water bath, the SV screen displays the command “dHAL”, and the ALM indicator flashes (red). To make the alarm stop, turn OFF the equipment.

To change the timer’s value, proceed the following way:

- Press SET for 3 seconds. The PV screen will display “Loc” and the SV screen will display, by default, the digits “0000”.
- Use ▲ and ▼ to configure the “0003” option and press SET. The PV screen will display the command “dHAL”.
- Use ▲ and ▼ (and ◀ to shift digits) to set the alarm’s temperature. Press SET to save and wait to go back to the initial screen.

## 6.4. Draining

All the water baths, except TWBA-001-001 series, have a draining system, which includes a pump and a drainpipe.

- Before draining, make sure the ON/OFF switch is in the OFF position (O); otherwise, the resistor could be functioning once the bath has been emptied and get damaged.
- With the equipment plugged, press the drain switch; the drainpipe will start emptying the water bath container.

### 6.5. Autosetting

The autosetting function allows you to adjust the internal parameters at a certain working temperature. To do that:

- Select the working temperature (follow instructions in 6.1).
- Press ▼ until the AT indicator starts flashing (Green) –you just selected the autosetting function.
- The autosetting will end once the AT indicator turns off.
- After the autosetting, the equipment will optimize its service.

### 6.6. “Stop” Function

The “Stop” function allows you to stop and resume the equipment throughout the execution of a function. Through this function, you can stop and resume the equipment while heating and if the timer is on, you can stop and resume it.

- In the initial screen, press ▲ for a few seconds. You will hear a “click” and the SV screen will display “Stop”.
- To resume the function, press ▲ for a few seconds again until you hear a light beep that indicates the resumption.

## 7. Fault Resolution

- The equipment cannot be turned ON
  - Check whether the power cable is plugged in.
  - Make sure the cable is rightly connected.
- The temperature cannot reach the set point
  - Check for apparent damages that may have arisen during transportation
  - The resistor and/or the electric system may be damaged.

*If these faults are not resolved, contact your manufacturer/supplier.*

## 8. Cleaning and maintenance

- First, disconnect the equipment from the socket during cleaning and/or maintenance.
- Proper maintenance can keep instruments working and lengthen their lifetime.
- Do not spray cleanser onto the instrument when cleaning.
- The equipment must be cleaned and disinfected when before shipping, and packed into its original packaging.
- Always use the equipment in a dry, clean space with a stable temperature.

## 9. Storage and Transportation

- Keep the equipment in a dry, clean place with good with good ventilation and no corrosive gases.
- During transportation, prevent the equipment from wetting and avoid violent collision.

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