

**GENERAL INFORMATION**
**Product name : Scissors**
**Description :** Made of stainless steel. Not suitable for autoclaving


**TECHNICAL DATA**

reference	l(mm)	type	pcs/pack
SCIS-14A-001	145	dissection straight, pointed/pointed (A)	1
SCIS-14B-001	145	dissection straight, pointed/blunt (B)	1
SCIS-14C-001	145	dissection straight, blunt/blunt (C)	1
SCIS-115-001	115	microscopy straight, pointed/pointed (D)	1
SCIS-220-001	220	for general use, with PVC handle (E)	1

**PACKAGING AND LOGISTICS**

reference	vol (l)	kg	TARIC	GTIN
SCIS-14A-001	0,056	0,046	90278099	08434868041288
SCIS-14B-001	0,048	0,048	90278099	08434868041301
SCIS-14C-001	0,048	0,048	90278099	08434868041325
SCIS-115-001	0,042	0,022	90278099	08434868041264
SCIS-220-001	0,384	0,88	90278099	08434868041349

PRODUCT PHOTO



## INFORMACIÓN GENERAL

**Nombre del producto :** Tijeras**Descripción :** De acero inoxidable. No apto para autoclave


## DATOS TÉCNICOS

referencia	l (mm)	tipo	unidades por ref.
SCIS-14A-001	145	disección recta, aguda/aguda (A)	1
SCIS-14B-001	145	disección recta, aguda/roma (B)	1
SCIS-14C-001	145	disección recta, roma/roma (C)	1
SCIS-115-001	115	microscopía recta, aguda/aguda (D)	1
SCIS-220-001	220	uso general, mango de PVC (E)	1

## EMBALAJE Y DATOS LOGÍSTICOS

referencia	vol (l)	kg	TARIC	GTIN
SCIS-14A-001	0,056	0,046	90278099	08434868041288
SCIS-14B-001	0,048	0,048	90278099	08434868041301
SCIS-14C-001	0,048	0,048	90278099	08434868041325
SCIS-115-001	0,042	0,022	90278099	08434868041264
SCIS-220-001	0,384	0,88	90278099	08434868041349

FOTO DEL PRODUCTO



## INFORMATIONS GÉNÉRALES


**Nom produit : Ciseau****Description :** En acier inoxydable. Ne convient pas à l'autoclavage

## DONNÉES TECHNIQUES

référence	long. (mm)	type	unités par ref.
SCIS-14A-001	145	à dissection droit, pointu/pointu (A)	1
SCIS-14B-001	145	à dissection droit, pointu/arrondi (B)	1
SCIS-14C-001	145	à dissection droit, arrondi/arrondi (C)	1
SCIS-115-001	115	à microscopie droit, pointu/pointu (D)	1
SCIS-220-001	220	à usage général, manche en PVC (E)	1

## EMBALLAGE ET LOGISTIQUE

référence	vol (l)	kg	TARIC	GTIN
SCIS-14A-001	0,056	0,046	90278099	08434868041288
SCIS-14B-001	0,048	0,048	90278099	08434868041301
SCIS-14C-001	0,048	0,048	90278099	08434868041325
SCIS-115-001	0,042	0,022	90278099	08434868041264
SCIS-220-001	0,384	0,88	90278099	08434868041349



## INFORMAZIONE GENERALE

**Nome del prodotto :** Forbici**Descrizione :** In acciaio inox. Non adatto per la sterilizzazione in autoclave

## DATI TECNICI

referenza	L (mm)	tipo	unità per ref.
SCIS-14A-001	145	per dissezione rette, acuta/acuta (A)	1
SCIS-14B-001	145	per dissezione rette, acuta/smussa (B)	1
SCIS-14C-001	145	per dissezione rette, smussa/smussa (C)	1
SCIS-115-001	115	per microscopia rette, acuta/acuta (D)	1
SCIS-220-001	220	per uso generale, manico in PVC (E)	1

## IMBALLAGGIO E DATI LOGISTICI

referenza	vol (l)	kg	TARIC	GTIN
SCIS-14A-001	0,056	0,046	90278099	08434868041288
SCIS-14B-001	0,048	0,048	90278099	08434868041301
SCIS-14C-001	0,048	0,048	90278099	08434868041325
SCIS-115-001	0,042	0,022	90278099	08434868041264
SCIS-220-001	0,384	0,88	90278099	08434868041349

FOTO DEL PRODOTTO

